


DxC600/CX INSTRUMENTS CALIBRATION VERIFICATION - 9/2006

Synchron Control Level 1 lot #M511371	Vigil Serology Control Level 1 lot#M504291	D1 Control lot#M412331 MHC/Fin
Synchron Control Level 2 lot #M511372	Vigil Serology Control Level 2 lot#M409332	D1 Control lot#M505211 Dyersville
Synchron Control Level 3 lot #M511373	Vigil Serology Control Level 3 lot#M504293	D2 Control lot#M505212
Amm/Alc Control Level 1 lot #M508421 Fin		D3 Control lot#M505213
Amm/Alc Control Level 1 lot #M508421 MHC		
Amm/Alc Control Level 2 lot #M508422 Fin		
Amm/Alc Control Level 2 lot #M508423 MHC		
Amm/Alc Control Level 3 lot #M406493	HbA1c Calibrator 2 lot #M602005	T4 Low Cal lot#56384229
Biorad Multiquel Liquid Assayed Level 1 lot#45531		T4 High Cal lot #56384229
Biorad Multiquel Liquid Assayed Level 2 lot#45532		
Biorad Multiquel Liquid Assayed Level 3 lot#45533		

Again, the tolerances are pretty wide but, with the exception of the ACL 9000 Protime, recovered results are running to the middle of the ranges. More importantly, we continue to show excellent intersite agreement.  
 sraymond 20061010



Ranges were obtained from control package inserts for the Beckman DxC Instruments.

Beckman Synchron Controls																			
Site	MAE	MAW	MHC	Finley	Dyer		MAE	MAW	MHC	Finley	Dyer		MAE	MAW	MHC	Finley	Dyer		
Level	Level 1	Level 1	Level 1	Level 1	Level 1	Level 1	Level 2	Level 2	Level 2	Level 2	Level 2	Level 2	Level 3	Level 3	Level 3	Level 3	Level 3	Level 3	Level 3
Assay	Result	Result	Result	Result	Result	Range	Result	Result	Result	Result	Result	Range	Result	Result	Result	Result	Result	Result	Range
Albumin	NA	NA	2.3	2.3	2.3	1.9-2.7	NA	NA	3.8	3.7	3.7	3.3-4.1	NA	NA	5.2	5.2	5.2	5.2	4.5-5.7
ALP	NA	NA	40	41	39	30-50	NA	NA	155	155	151	128-178	NA	NA	270	267	261	226-306	
ALT	NA	NA	16	16	17	13-23	NA	NA	163	163	157	137-197	NA	NA	308	306	293	284-348	
Amylase	NA	NA	45	43	44	36-56	NA	NA	219	220	220	193-253	NA	NA	394	395	391	362-442	
AST	NA	NA	21	22	22	17-27	NA	NA	197	196	189	173-223	NA	NA	354	350	340	300-420	
BUN3	7	7	NA	NA	NA	5-9.0	33	33	NA	NA	NA	28-36	60	61	NA	NA	NA	52-64	
BUN	NA	NA	8	8	8	6.0-10.0	NA	NA	36	36	36	32-40	NA	NA	64	62	64	57-69	
Calcium	7.8	8.1	7.6	7.6	7.5	7.1-8.3	10.6	10.9	10.6	10.3	10.5	9.8-11.4	13.3	13.5	13.2	13	13.3	12.1-14.5	
Chloride	86	83	80	81	81	77-85	104	103	101	101	101	97-107	122	122	118	120	118	115-127	
Cholesterol	NA	NA	100	103	98	89-119	NA	NA	163	165	167	148-168	NA	NA	229	228	216	210-250	
CK	NA	NA	56	56	55	48-68	NA	NA	351	369	361	327-427	NA	NA	637	661	642	579-779	
CO2	11	11	12	12	11	9.0-15.0	21	20	21	21	20	17-23	30	29	29	30	29	24-34	
CRE3	0.5	0.5	NA	NA	NA	0.3-0.7	4	3.9	NA	NA	NA	3.7-4.5	7.5	7.3	NA	NA	NA	6.9-8.5	
Creatinine	NA	NA	0.5	0.5	0.6	0.4-0.8	NA	NA	3.9	4	4	3.6-4.4	NA	NA	7.1	7.3	7.3	6.4-8.0	
Iron	NA	NA	NA	51	NA	37-63	NA	NA	NA	171	NA	151-191	NA	NA	NA	279	NA	248-320	
GGT	NA	NA	14	15	16	7.0-27.0	NA	NA	1	173	162	137-197	NA	NA	326	326	306	259-379	
Glucose	NA	NA	NA	NA	45	37-53	NA	NA	224	NA	232	216-256	NA	NA	NA	NA	412	384-452	
GLUm	42	40	41	44	NA	32-56	226	218	NA	222	NA	215-241	385	370	395	388	NA	365-435	
HDL	NA	NA	27	27.4	26	21-33	NA	NA	47.8	47.2	47.1	38-52	NA	NA	68.5	66.5	69.1	55-73	
IBCT	NA	NA	NA	216	NA	150-230	NA	NA	NA	359	NA	262-382	NA	NA	NA	494	NA	360-516	
Potassium	2.5	2.5	2.5	2.4	2.5	2.2-2.8	5.1	5.1	5	4.9	5.1	4.6-5.4	7.7	7.7	7.5	7.5	7.4	6.9-8.1	
LD-L	NA	NA	44	49	NA	33-63	NA	NA	185	196	NA	168-248	NA	NA	319	335	NA	299-419	
LDLD	NA	NA	42.8	43.5	43.2	36-48	NA	NA	76.4	74.2	74.2	65-83	NA	NA	110.6	100	111.3	92-126	
Lipase	NA	NA	120	124	NA	96-136	NA	NA	72	76	NA	53-85	NA	NA	25	24	NA	14-34	
Lithium	NA	NA	0.56	NA	NA	0.4-0.8	NA	NA	1.51	NA	NA	1.2-1.8	NA	NA	2.42	NA	NA	1.9-2.9	
Magnesium	NA	NA	1.22	1.18	NA	0.8-1.6	NA	NA	2.31	2.26	NA	2.0-2.8	NA	NA	3.53	3.44	NA	2.9-4.1	
Salicylate	NA	NA	30.5	NA	NA	23-37	NA	NA	18.3	NA	NA	14-24	NA	NA	8.8	NA	NA	5-11.0	
Sodium	114	117	115	115	115	110-120	141	143	140	139	141	134-146	167	167	163	165	164	157-171	
Phosphorus	NA	NA	2	2	NA	1.6-2.4	NA	NA	4.5	4.3	NA	3.7-5.1	NA	NA	7.1	6.6	NA	5.8-7.8	
Total Bilirubin	NA	NA	0.9	1	0.9	0.5-1.5	NA	NA	4.3	4.2	4.1	3.0-5.2	NA	NA	7.4	7.5	7.7	5.4-9.4	
Total Protein	NA	NA	3.8	4	3.9	3.6-4.4	NA	NA	6.1	6.1	6.1	5.6-6.8	NA	NA	8.3	8.4	8.2	7.4-9.2	
Triglyceride	NA	NA	67	68	70	56-83	NA	NA	107	107	106	90-130	NA	NA	148	141	143	123-173	
Uric Acid	NA	NA	2.5	2.5	NA	2.2-3.0	NA	NA	6.9	6.8	NA	5.9-7.5	NA	NA	11.1	11	NA	9.6-11.6	

DxC600/CX INSTRUMENTS CALIBRATION VERIFICATION - 9/2006

Beckman Ammonia/Alcohol Controls									
Site	MHC	Finley		MHC	Finley		MHC	Finley	
Level	Level 1	Level 1	Level 1	Level 2	Level 2	Level 2	Level 3	Level 3	Level 3
Assay	Result	Result	Range	Result	Result	Range	Result	Result	Range
Ammonia	85	89	48-116	254	250	197-265	1107	1153	1040-1193
Alcohol	0.05	0.05	0.04-0.06	0.11	0.11	0.09-0.12	0.49	0.48	0.421-0.509

Excellent

Beckman Vigil Serology Controls						
Site	Finley		Finley		Finley	
Level	Level 1	Level 1	Level 2	Level 2	Level 3	Level 3
Assay	Result	Range	Result	Range	Result	Range
ASO	87.1	62.3-101.5	200.2	166-226	383.2	332-422

Beckman Vigil Serology Controls						
Site	MHC		MHC		MHC	
Level	Level 1	Level 1	Level 2	Level 2	Level 3	Level 3
Assay	Result	Range	Result	Range	Result	Range
CRP	1	0.7-1.7	4.6	4.5-5.7	7.2	6.5-8.5

Beckman T4 Calibrators						
Site	MHC					
Level	Low Cal	Low Cal	33%	67%	High Cal	High Cal
Assay	Result	Target	Result	Result	Result	Target
T4	0.01/0.0	0.2	6.1/5.7	13.0/13.3	20.2/19.6	19.7

Falls off a bit as [ ] decreases. 13.2, 6.5. Not typical & not significant. Will wait on next set of challenges.

Beckman D1/D2/D3 Bilirubin Controls												
Site	MHC	Finley	Dyer		MHC	Finley	Dyer		MHC	Finley	Dyer	
Level	D1	D1	D1	D1	D2	D2	D2	D2	D3	D3	D3	D3
Assay	Result	Result	Result	Range	Result	Result	Result	Range	Result	Result	Result	Range
Direct Bilirubin	0.9	0.9		0.5-1.3 MHC/F	3.9	3.8	3.8	3.4-4.4	5.7	5.6	5.7	5.1-6.5
			0.7	0.4-1.2 Dyers								

Beckman Synchron Level 1 Control diluted 1:10 (Dilution prepared using Saline. Run samples as a CSF sample.)										
Site	MHC	Finley	MHC	Finley	MHC	Finley	MHC	Finley	MHC	Finley
Dilution	Saline	Saline	30%	30%	50%	50%	70%	70%	100%	100%
MTP	RSL	RSL	95	101	162	172	222	223	320	RSH

OK. Would like to look for some MTP-assayed material. ...Lynn

Biorad Multiquel Liquid Assayed Controls (Acid Phos slope set to 1.0)						
Site	MHC					
Assay	Level 1	Level 1 Range	Level 2	Level 2 Range	Level 3	Level 3 Range
Acid Phos	5.6	5.5-8.3	14.2	13.1-19.7	31.3	26.9-40.3

Good.

Beckman HbA1c Cal 2			
Site	Dyer		
Level	Hem Rgt	1:2 Dil	HbA1c Cal 2
Assay	Result	Result	Target
Hb	RSL	8.2	17.6
			17.8

BNP Calibration Verification kit lot #162422	MHC BNP Cartridge lot #W38825B	I-Stat Cardiac Marker L1 lot #M061591
i-Stat Calibration Verification Material lot #B0656830	Fin BNP Cartridge lot #W38623B	I-Stat Cardiac Marker L2 lot #M061592
MAE i-Stat cartridge lot #Y05266	Beckman D1 Bilirubin Control lot #M412331	I-Stat Cardiac Marker L3 lot #M061593
MHC i-Stat cartridge lot #606095	Beckman D2 Bilirubin Control lot #M505212	MHC I-Stat1 Troponin I cartridge lot #U06131
Dyer i-Stat cartridge lot #U06095/J06117	Beckman D3 Bilirubin Control lot #M505213	Finley I-Stat1 Troponin I cartridge lot#U0617A
Dyer BNP cartridge lot#W38720B	Beckman D4 Bilirubin Control lot #M505214	Wescor Sweat Controls lot#1504237

Site	MHC	Finley	Dyers		MHC	Finley	Dyers		MHC	Finley	Dyers		MHC	Finley	Dyers		MHC	Finley	Dyers	
Level	Level A	Level A	Level A	Level A	Level B	Level B	Level B	Level B	Level C	Level C	Level C	Level C	Level D	Level D	Level D	Level D	Level E	Level E	Level E	Level E
Assay	Result	Result	Result	Range	Result	Result	Result	Range	Result	Result	Result	Range	Result	Result	Result	Range	Result	Result	Result	Range
BNP	20.9	19	22.2	<5.0-33.6	39.5	34	46.4	17.9-60.1	945	709	991	492-1270	3830	3560	4620	1750->5000	>5000	3420	>5000	1990->5000

Site	MHC	Dyer	MAE		MHC	Dyer	MAE		MHC	Dyer	MAE		MHC	Dyer	MAE		MHC	Dyer	MAE	
Level	Level 1	Level 1	Level 1	Level 1	Level 2	Level 2	Level 2	Level 2	Level 3	Level 3	Level 3	Level 3	Level 4	Level 4	Level 4	Level 4	Level 5	Level 5	Level 5	Level 5
Assay	Result	Result	Result	Range	Result	Result	Result	Range	Result	Result	Result	Range	Result	Result	Result	Range	Result	Result	Result	Range
pO2	67	67	67	50.5-76.5	76	82	85	62.8-88.8	96	101	100	82.3-106.9	123	123	120	100.9-130.9	353	343	362	296.7-400.3
pCO2	89	89.6	85.5	79.4-96.64	62	61.8	60.8	54.86-67.10	26	26.1	25.3	20.40-31.18	18	17.6	16.7	11.87-22.15	16	15.5	15	10.47-20.27
pH	6.825	6.834	6.826	6.779-6.877	7.171	7.18	7.163	7.126-7.224	7.45	7.456	7.457	7.409-7.507	7.697	7.706	7.706	7.660-7.758	8.112	8.141	8.147	8.095-8.193

Site	MAE	MAW		MAE	MAW		MAE	MAW		MAE	MAW	
Level	Level 1	Level 1	Level 1	Level 2	Level 2	Level 2	Level 3	Level 3	Level 3	Level 4	Level 4	Level 4
Assay	Result	Result	Range	Result	Result	Range	Result	Result	Range	Result	Result	Range
Total Bilirubin	1.6	1.7	1.2-2.0	5.9	6.2	5.2-6.6	8.6	8.9	7.9-9.7	16.6	17	no ranges
Direct Bilirubin	0.7	0.5	0.5-1.3	3.7	3.8	3.4-4.4	5.5	5.3	5.1-6.5	11	11.2	no ranges

Site	MHC					
Level	Level 1	Level 1	Level 2	Level 2	Level 3	Level 3
Assay	Result	Range	Result	Range	Result	Range
Chloride	43	38-46	73	65-79	137	121-147

Site	MHC	Finley		MHC	Finley		MHC	Finley	
Level	Level 1	Level 1	Level 1	Level 2	Level 2	Level 2	Level 3	Level 3	Level 3
Assay	Result	Result	Range	Result	Result	Range	Result	Result	Range
Troponin I	0.25	0.26	0.20-0.36	1.89	1.6	1.27-2.37	34.23	36.26	20.61-54.35

September 2006 LM

Within tolerances but the tolerances are very forgiving. Finley's cartridge lot is generating lower values relative to the others. May consider increased management of cartridge lots. With a neg predictive value cutoff of <100 pg/ml and the fact that the bnp is run "within house", the difference expressed here at very high levels would have no clinical impact. Will take a look at logistics of cartridge lot management anyway.

Hemosil Assayed Low Abnormal Control lot #N0957718	MHC D-Dimer Low Control lot #B60418
Hemosil Assayed Normal Control lot #N1050340	MHC D-Dimer High Control lot #B60418
Hemosil Assayed High Abnormal Control lot #N0160794	Finley D-Dimer Low Control lot #B50363
Hemosil Assayed Low Fibrinogen Control lot #N0261075 MHC	
Hemosil Assayed Low Fibrinogen Control lot #N1150363 Fin	Finley D-Dimer High Control lot #B50363

ACL 9000								
Hemosil Assayed Control Material								
Site	MHC							
Level	Low Abn	Low Abn	Normal	Normal	High Abn	High Abn	Low Fib	Low Fib
Assay	Result	Target	Result	Target	Result	Target	Result	Target
Prottime	36.4	27.2-36.8	11.2	9.1-12.1	60.7	42.3-57.3	NA	NA
Fibrinogen	194	142-212	356	272-392	NA	NA	121	87-125
Thrombin	15.4	14.5-19.7	12.2	12.2-16.2	NA	NA	NA	NA

Low Abn & High Abn a little high to the target on this set. Wasn't the case last time. Repeat and we'll take it from there.

ACL Classic											*ACL 9000 range							
Hemosil Assayed Control Material											no range listed for ACL Classic							
Site	Finley	Dyer	MAE	MAW		Finley	Dyer	MAE	MAW		Finley	Dyer	MAE	MAW		Finley		
Level	Low Abn	Low Abn	Low Abn	Low Abn	Low Abn	Normal	Normal	Normal	Normal	Normal	High Abn	High Abn	High Abn	High Abn	High Abn	Low Fib	Low Fib	
Assay	Result	Result	Result	Result	Target	Result	Result	Result	Result	Target	Result	Result	Result	Result	Target	Result	Target	
Prottime	29.6	31.5	30.3	30.9	24.1-32.5	11	10.9	11.5	10.3	8.8-11.8	51.6	52.5	52.6	52.3	42.3-57.3	NA	NA	
Fibrinogen	186	NA	NA	NA	139-179	363	NA	NA	NA	287-377	NA	NA	NA	NA	NA	124	81-131	

Very good

ACL 9000						
D-Dimer						
Site	MHC					
Level	Low Control	Low Control	1:1 dil	Calc target	High Control	High Control
Assay	Result	Range	Result	Range	Result	Range
D-Dimer	306	267-445	519	516	739	575-777

Good!

ACL 7000						
D-Dimer						
Site	Finley					
Level	Low Control	Low Control	1:1 dil	Calc target	High Control	High Control
Assay	Result	Range	Result	Range	Result	Range
D-Dimer	338	247-411	483	510	665	587-795